

NARA® LURE

Allergen-free monitor for rats and mice

Allergen-management has become a major topic in many food product companies.

NARA® is the first allergen free monitor!

Law as well as many food standards now prohibit the use of anticoagulants for monitoring purposes. Only after the determination of rodent activity anticoagulants may be used subject to stringent restrictions. Common monitors have the critical disadvantage that they serve as food for rodents, which increases their lifetime and the risk of a rising reproduction of rodents. Moreover, most common baits are prone to mould formation in damp areas. In other areas they are often infested by non-target organisms, e.g. flour beetles, cockroaches, bristletails, moths, ants or snails.

No other beats NARA® Lure: 3 up to 12 months long lasting efficacy in 3 different aromas will be your professional aid for long-term monitoring. Rely upon award winning NARA® quality for guaranteed success in any environment.

With the new aromas there is no other product that can deliver this level of quality and flexibility. **Choose between vanilla, meat or fish aroma.**

Tuna or peanut butter leave a big mess that has to be cleaned up after some days of usage. Professionals now rely upon NARA® Lure, which is the most hygienic and long-lasting way to monitor. No cleaning up afterwards and no fast decrease of efficient monitoring with up to 12 months of non-stop monitoring. Only with NARA®!

Packaging unit 100 pieces/jar



Vanilla flavour

Meat/Beef flavor

Fish flavour

Just place the Lure into the bait cup of the unset Snap Trap (use high quality Gorilla Traps, Kness or similar) or position in any station if you want to monitor activity before it can cause damage.

Remember: monitoring will only be successful if all the pre-work is done: proofing, sealing, hygiene measurements, funneling, using the correct stations (e.g. RUNBOX)... Ask us for support!